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## Web portal ID

<http://easimage.imadev.ch>

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A screenshot of the login page for the LACOLLINE RADIOLOGY PORTAL. The page has a dark blue header with the logo of 'IMAGERIE & DEVELOPPEMENT' and the text 'LACOLLINE RADIOLOGY PORTAL'. Below the header, there is a login form with two input fields: 'Username:' and 'Password:'. Below the password field is a 'Login' button and a link for 'Password Forgotten?'. The page is otherwise blank.

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The screenshot shows the main interface of the LACOLLINE RADIOLOGY PORTAL. At the top, there is a dark blue header with the logo on the left, the text "LACOLLINE RADIOLOGY PORTAL" in the center, and the user's login status "Logged in: [1] narc | Logout" on the right. The number "1" is enclosed in a box. Below the header, the "SEARCH" section contains three input fields: "Name:", "Patient ID:", and "Accession Number:", each with a "Search" button. The number "2" is enclosed in a box next to these fields. Below the search section, the "BROWSE" section lists various filters: "Last 6 hours", "Today", "Study List", "Just Added", "Last Week", "Last Week CT", "Last Week MR", "Today CT", "Today MR", "Yesterday", "Yesterday CT", and "Yesterday MR". The number "3" is enclosed in a box next to the "Last Week" filter. At the bottom of the page, a small footer reads "Service provided by ID Imagerie & Développement."

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The screenshot shows a web interface for medical imaging search results. At the top, there is a header with the logo 'IMAGERIE & DEVELOPPEMENT' and 'SEARCH RESULTS'. The user is logged in as 'admin' with a 'Logout' button. Below the header, there is a navigation link '< HOME' and a section '3 Studies Found'. The main content area displays a list of studies with the following details:

Study Name	Date/Time	Actions
AGECANONIX Specials 1Co... s, OQYOrzL, 0 with_spiral_CTA_pre	12.12.03 09:33	  
BRAINIX MR, 7 series, 5Yp0E, 0 IRM cérébrale, neuro-crâne	01.12.06 14:39	  
KNIX MR, 6 series, ozp005jY2xG Knee (R)	01.01.07 12:00	  

Numbered callouts on the screenshot indicate the following actions:

- 1: Clicking on the study name.
- 2: Clicking on the OsiriX icon.
- 3: Clicking on the Download icon.
- 4: Clicking on the Weasis icon.

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
IMAGERIE & DEVELOPPEMENT  
BRAINIX - IRM CÉRÉBRALE, NEURO-CRÂNE  
Logged in as admin | Logout

SEARCH RESULTS FOR: NIX

**STUDY INFORMATION**

Patient's name: BRAINIX  
Date of birth: 01.03.49  
Patient ID: 5Yp0E  
Accession number: 0  
Study date: 01.12.06 14:39  
Modality: MR  
Description: IRM cérébrale, neuro-crâne

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
---

Show All Frames

Viewing this series requires Adobe Flash Player

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